

# ***SBEAMS-Microarray*** **User Guide**

Bruz Marzolf  
April 12, 2006


## Table Of Contents

<b>1</b>	<b><i>Basic Data Access</i></b>	<b>3</b>
1.1	Getting Into <i>SBEAMS</i>	3
1.2	Download Array Data	3
1.2.1	Access the Download Page	3
1.2.2	File Types	4
1.3	Annotate Affy Samples or Arrays	5
1.3.1	Annotate Sample Page	5
1.3.2	Annotate Array Page	6
1.4	Simple Query	6
1.5	Advanced Query	7
1.5.1	Entering Query	7
1.5.2	Viewing Query Results	7
<b>2</b>	<b><i>Analysis Pipeline</i></b>	<b>7</b>
2.1	Normalization	7
2.1.1	Starting a New Analysis Session	7
2.1.2	Selecting and Grouping Arrays	8
2.1.3	Choosing Normalization Options	10
2.1.4	Viewing Results	11
2.2	Differential Expression Testing	11
2.2.1	Return to Normalization Set	11
2.2.2	Perform SAM Run	11
2.2.3	Adding data to GetExpression	12
2.3	Querying GetExpression and Launching Cytoscape	13
2.3.1	Query Data	13
2.3.2	Launch Cytoscape	13
2.3.3	Viewing Expression Data in Cytoscape	14

# 1 Basic Data Access

## 1.1 Getting Into *SBEAMS*

- 1) Login:
  - a. <http://<SBEAMS server name>/sbeams>



**Figure 1.1** Login page for *SBEAMS*

- b. Click on the link **Login to SBEAMS**
    - c. Login with your username and password as set up by your *SBEAMS* Administrator
  - 2) Go to Microarray Module
    - d. On the left menu bar click on **Microarray** to enter *SBEAMS-Microarray*

## 1.2 Download Array Data

Many array data files can be downloaded or viewed directly from *SBEAMS*

### 1.2.1 Access the Download Page

To download or view array files from a particular project choose the **Download Data** link in the left menu.




Figure 1.2 Download Data page

- 1) Select a project from the drop down list located at the top of the page.
- 2) To download files, select the check boxes for the file type you'd like for each array you'd like to download. To download one type of file for all arrays in the project, check the box next to the file type at the top of the table (e.g. CEL, RPT). The default setting is to download all the CEL files from a project.
- 3) Once all the files are selected click the button **GET\_AFFY\_ARRAY\_FILES**
- 4) Select a location to save the zip file containing all the information selected for downloading

## 1.2.2 File Types

Here is a brief description of all the file types available for viewing or downloading from SBEAMS

File Extension	Description	View File Directly	Download
CHP	CHP. Binary Affymetrix file. CHP files contain probe set analysis results generated from Affymetrix software.	NO	Yes
CEL	CEL. Binary Affymetrix file. The CEL file stores the results of the intensity calculations on the pixel values of the DAT file	NO	Yes
XML	XML. MAGE XML Affymetrix file. Contains information from Affymetrix GCOS Software collected during sample preparation, hybridization, washing and scanning.	Yes	Yes

RPT	RPT. Text report. Contains information about the CHP file, used for basic quality control	Yes	Yes
R_CHP	R_CHP. Text File. Contains Probe set intensity values, calculated by using R/Bioconductor affy mas5.0 algorithms	No	Yes
JPEG	JPEG. Jpeg image of the Affy Chip generated by R using the image method within the affy library	Yes	Yes
EGRAM_PF.jpg	EGRAM_PF.jpg. Electrophoregram image of the Pre-fragmented cRNA	Yes	Yes
EGRAM_T.jpg	EGRAM_T.jpg. Electrophoregram image of the total RNA	Yes	Yes
EGRAM_F.jpg	EGRAM_F.jpg. Electrophoregram image of the fragmented cRNA	Yes	Yes

### 1.3 Annotate Affy Samples or Arrays

All samples hybridized to Affymetrix arrays can be annotated and viewed within SBEAMS.

#### 1.3.1 Annotate Sample Page

To access the page to annotate or view information about a sample, choose the **Microarray Home** link in the left menu.

The screenshot shows the SBEAMS - MicroArray web interface. The main content area is titled "Maintain Affy Array Sample". At the top, there is a login section with "Login: bmarzoll (3)" and "Group: Microarray\_admin". Below this, there is a "Project" dropdown menu set to "bmarzoll - ISB Test Data (2)". The form fields are as follows:

- Project:** bmarzoll - ISB Test Data (2)
- Sample Tag:** LN\_1
- Full Sample Name:** (empty)
- Sample Group Name:** LN
- Sample Provider:** Institute for Systems Biology
- Organism:** Human
- Strain or Line:** (empty)
- Individual:** (empty)
- Sex:** (empty)
- Age:** (empty)

Figure 1.3 Annotating a Sample

- 1) Select a project from the drop down located at the top of the page.
- 2) Click on a Link under the "Sample\_Tag" Column to view information about a sample

- 3) Enter Data to describe the sample
- 4) To save time entering multiple samples, use the "Save Template" functionality. This will save a copy of the current page. To re-use the template on your next sample, choose the template in the "Existing Templates" drop down box. Then click **Set Fields to this Template**. Change any fields that are different for the second sample.
- 5) Click "Insert" if this is a new record or "Update" if the record was modified to save the information

### 1.3.2 Annotate Array Page

To access the page to annotate or view information about a specific Affy Array, choose the **Microarray Home** link on the left menu.

- 1) Select the project from the drop down located at the top of the page.
- 2) Click on a Link under the "File\_root" Column to view information about an array
- 3) Enter Data to describe the array
- 4) Click "Insert" or "Update" to save the information

### 1.4 Simple Query

Perform simple queries within a project to view probe set intensity values. Expression values are taken from the R\_CHP files.

Choose the **Affy Gene Intensity** link on the left menu to access the Simple Query page.




Figure 1.4 Simple Query interface

Enter a query with an Affy Probe Set ID, Gene Symbol, Gene Name or Accession number:

- 1) Choose the project of interest
- 2) Select specific arrays, default is to have all arrays for the project selected
- 3) Enter the search term
- 4) Hit the "simply\_query" button

## 1.5 Advanced Query

Perform Advanced queries within one or many projects to view probe set intensity values. Expression values are taken from the R\_CHP files.

### 1.5.1 Entering Query

- 1) Change to the Advanced Query page, select one or more projects, arrays and query constraints:
- 2) Choose the project(s) of interest
- 3) Select the arrays in the chosen projects you wish to query
- 4) Enter search constraints such as the probe set ID, gene symbol
- 5) Select the Data Columns you want to display, such as Signal Intensity and/or Detection P-value
- 6) Select Display Options
- 7) Hit the "QUERY" button to run your query

### 1.5.2 Viewing Query Results

- Data can be displayed in a html table, tsv, csv, excel or xml formats
- Any of the columns may be sorted
- Link to Affy annotation page is provided

## 2 Analysis Pipeline

### 2.1 Normalization

Using the pipeline begins with starting a new analysis session, choosing and grouping arrays, and then picking normalization options.

#### 2.1.1 Starting a New Analysis Session

- 1) Switch to a project of interest using the drop-down list at the top of the page.
- 2) Choose the **Data Pipeline** link on the left menu bar
- 3) Start a new Analysis Session
  - a. Click "Start Session"

## 2.1.2 Selecting and Grouping Arrays

The screenshot shows the SBEAMS - MicroArray web interface. The top header displays the user login as 'bmarzolf (3)' with the group 'Microarray\_admin' and the project 'kstegmaier - Primary APL Samples (1)'. Below the header, there are tabs for 'File Groups', 'Normalized Data', and 'Analysis Results'. The 'File Groups' tab is active, and the 'Start a New Analysis Session' button is visible. The main content area is titled 'Select Additional Projects To view arrays to include in analysis' and contains a dropdown menu with two options: 'bmarzolf - ISB Test Data - (2)' and 'kstegmaier - Primary APL Samples - (1)'. Below this, there is a section titled 'Please Select the arrays to utilize in the analysis pipeline' with a 'Click to select or de-select all arrays' instruction. A table of arrays is displayed with columns for File\_Root, Affy Chip Design, Sample\_Tag, Sample Group Name, Full\_Sample\_Name, Organism, and CEL. The table contains 15 rows of data, each with a checkbox in the CEL column. At the bottom of the table, there are 'Add Arrays' and 'Reset' buttons.

File_Root	Affy Chip Design	Sample_Tag	Sample Group Name	Full_Sample_Name	Organism	CEL
20030609_001_DP2002050118AA	HG-U133A	AML2	primary_patient_AML_cells		Human	<input checked="" type="checkbox"/>
20030609_001_CL2002042638AA	HG-U133A	NEU2	normal_human_neutrophils		Human	<input checked="" type="checkbox"/>
20030609_001_DP20020612PK7AA	HG-U133A	AML4	primary_patient_AML_cells		Human	<input checked="" type="checkbox"/>
20030609_001_CL2002042642AA	HG-U133A	MON3	normal_human_monocytes		Human	<input checked="" type="checkbox"/>
20030609_001_CL2002042640AA	HG-U133A	MON1	normal_human_monocytes		Human	<input checked="" type="checkbox"/>
20030609_001_CL2002042637AA	HG-U133A	NEU1	normal_human_neutrophils		Human	<input checked="" type="checkbox"/>
20030609_001_DP2002050119AA	HG-U133A	AML3	primary_patient_AML_cells		Human	<input checked="" type="checkbox"/>
20030609_001_DP20020612PK9AA	HG-U133A	AML5	primary_patient_AML_cells		Human	<input checked="" type="checkbox"/>
20030609_001_DP2002050117AA	HG-U133A	AML1	primary_patient_AML_cells		Human	<input checked="" type="checkbox"/>
20030609_001_CL2002042639AA	HG-U133A	NEU3	normal_human_neutrophils		Human	<input checked="" type="checkbox"/>
20030609_001_CL2002042641AA	HG-U133A	MON2	normal_human_monocytes		Human	<input checked="" type="checkbox"/>
20030609_001_DP20020612PK20AA	HG-U133A	AML6	primary_patient_AML_cells		Human	<input checked="" type="checkbox"/>

Figure 2.1 Adding files to a File Grouping

- 1) Multiple projects can be selected at the same time by clicking projects in the window "Select Additional Projects to view arrays to include in analysis"
- 2) Uncheck any arrays that should not be included in the analysis run.
- 3) Click the Add Arrays button at the bottom of the screen. This will present a list of the files you have selected.
- 4) If you wish to add more files to your list, click 'Show' at the bottom of the page to display the file adding interface. Add more arrays as described in steps 1-3.
- 5) Click the "Continue File Grouping" button



**SBEAMS - MicroArray**

db SBEAMS

SBEAMS Home  
Microarray Home  
Logout

● Login: **benarzell** (3) Group: Microarray\_admin  
● Project: kstegmaier-Primary APL Samples (1)

**Start a New Analysis Session**

File Groups | Normalized Data | Analysis Results

**Choose the number of Sample Comparison Groups**

**Sample Groups**

Group	Order	Sample Group Name	Reference Sample *
Sample Group 1	<input type="text" value="1"/>	<input type="text" value="normal_human_monocytes"/>	<input checked="" type="radio"/>
Sample Group 2	<input type="text" value="2"/>	<input type="text" value="normal_human_neutrophils"/>	<input type="radio"/>
Sample Group 3	<input type="text" value="3"/>	<input type="text" value="primary_patient_AML_cells"/>	<input type="radio"/>

The Reference Sample, will be compared to all additional samples groups provided if you wish to run t-test between two different sample groups. The "control group" should almost always be the Reference Sample, so that positive Log ratios indicate increased expression in the experimental group and vice versa.

Please Click "Update Order" if the Sample Group Names are changed

\* Please note that the reference sample can be ignored at the analysis so just two sample groups can be compared to one another.

SBEAMS - MicroArray 0.23-dev  
© 2005 Institute for Systems Biology

**Figure 2.2** Choosing a reference sample and sample group names in the File Grouping step

Add Sample Group information and select a reference sample:

- Replicate arrays will be combined together at the analysis phase if they are in the same sample group. The sample Group information is NOT utilized for any of the normalization calculations; it is simply collected here for use in further analyses that occur after normalization.
- 6) Change the number and names of groups as desired
  - 7) Click "Update Order"
  - 8) Associate CEL files with the appropriate Sample Group Names
  - 9) Click "Submit Group Names"
  - 10) Make sure all the files are in the correct sample groups
  - 11) Click "Complete File Grouping"

## 2.1.3 Choosing Normalization Options

The screenshot shows the SBEAMS MicroArray web interface. The main content area is titled "Affymetrix Expression Analysis: affy" and is in "Step 2: Select files for expression analysis:". Below this is a table with columns for "#", "File", and "Sample Name". The table contains 9 rows of data. Below the table, there are radio buttons for "RMA", "GC-RMA", and "Custom". The "Custom" option is selected. Under "Custom", there are dropdown menus for "Background Correction" (set to "rma"), "Normalization" (set to "quantiles"), "PM Correction" (set to "pmonly"), and "Summarization" (set to "medianpolish"). There are also checkboxes for "Log base 2 transform the results", "Produce MVA scatter plot", and "Produce correlation matrix". At the bottom, there are text input fields for "Enter description for analysis set (optional)" and "E-mail address where you would like your job status sent (optional)".

#	File	Sample Name
1	20030609_001_CL2002042640AA.CEL	MON1
2	20030609_001_CL2002042641AA.CEL	MON2
3	20030609_001_CL2002042642AA.CEL	MON3
4	20030609_001_CL2002042637AA.CEL	NEU1
5	20030609_001_CL2002042638AA.CEL	NEU2
6	20030609_001_CL2002042639AA.CEL	NEU3
7	20030609_001_DP2002050117AA.CEL	AML1
8	20030609_001_DP2002050118AA.CEL	AML2
9	20030609_001_DP2002050119AA.CEL	AML3

Figure 2.3 Normalization step options

- 1) You may edit the Sample names, which will be used to annotate the results
- 2) Select either RMA, GC-RMA or Custom
- 3) When Custom is chosen, you may separately choose the Background Correction, Normalization, PM Correction and Summarization options. These should be chosen with care as not all combinations are possible or desirable.
- 4) Optionally select whether Log base 2 results are desired (default, recommended)
- 5) Optionally select whether MVA plots and correlation matrices should be produced. These are useful as diagnostic plots. These require substantially more memory during analysis, and should be turned off to facilitate normalization of large numbers of arrays.
- 6) Optionally name the analysis (recommended to differentiate multiple analyses of the same data set)
- 7) Optionally enter an email address that will be sent a message upon completion of normalization run.
- 8) Click "Start Normalization"

## 2.1.4 Viewing Results

- 1) Either return to the browser window where you performed the analysis, or click on the result link in the email sent to you by the pipeline.
  - a. Click on the link **Show Files** located at the top of the page

The screenshot shows the SBEAMS - MicroArray web interface. At the top, there is a navigation menu with links for Home, Microarray Home, and Logout. Below this, there are login and project selection fields. The main content area is divided into three tabs: File Groups, Normalized Data, and Analysis Results. The Analysis Results tab is active, displaying 'Analysis Run Info' with fields for User Description, File Names, Analysis Description, and Start Additional Analysis. Below this, there is a table of 'Analysis Run Files' with columns for Show File, Download File, and Info. The table lists various files including 'Annotated expression values file', 'Tar Gzipped Archive of Analysis', 'Completed Job - Html File', '3ML file showing groupings', 'R Binary affy library expression file', 'R Error File', and 'R Script'.

Figure 2.4 Viewing available options and data in Normalization results

- 2) View the normalized data file
  - a. Under the header Analysis Run Files→Data click the link **Show**
  - b. This load log<sub>2</sub> expression values for all the genes and arrays just analyzed, although it's not very useful in the browser window.
- 3) To download the normalized data, under the header Analysis Run Files→Data click the link **Get**. These results may then be loaded into Excel or other software.

## 2.2 Differential Expression Testing

### 2.2.1 Return to Normalization Set

- 1) Make sure the project where you performed normalization is selected
- 2) Click **Data Pipeline** link on the left navigation menu
- 3) Click the "Normalization" tab
- 4) Click the **Show Files** Link for your normalization run

### 2.2.2 Perform SAM Run

- 1) Find the heading Start Additional Analysis->Multiple t-test
  - a. Click the link **Start Multitest**
- 2) Select Analysis method. Default is to run SAM

- a. Leave the radio button selection at “SAM”
- 3) Click the button “Next Step”
- 4) Review the reference sample and sample groups to be analyzed. Make sure the information is correct.

The screenshot shows the SBEAMS - MicroArray web interface. The main content area is titled "Multiple Testing: multtest" and "Step 3:". It contains several form fields and checkboxes. The "Reference Sample Group Name" is set to "C4-2". Below it, there is a checkbox labeled "Compare ONLY two sample group to one another" with a "YES" radio button. The "Information" section shows two sample groups: "Sample Group C4-2" with 2 files (20040421\_03\_C4-2\_2.CEL, 20040421\_04\_C4-2\_3.CEL) and "Sample Group LN" with 2 files (20040421\_01\_LN\_1.CEL, 20040421\_02\_LN\_2.CEL). The "Run SAM Analysis Two-class Unpaired Assuming Unequal Variances" section has a checked checkbox "Limit the HTML Results to FDR percent cut-off" set to 2%. Below this are two "AND" conditions: "A Minimum number of Genes" set to 10 and "A Maximum number of Genes" set to 250. There is also a checked checkbox "Include expression values in results". The "Name for analysis:" field is set to "LN vs. C4-2". An "E-mail address where you would like your job status sent. (optional)" field is empty. A "Submit Job" button is at the bottom.

Figure 2.5 Options for SAM Analysis

- a. When the run starts the reference sample will be compared to the other experimental sample groups.
- b. A single pairwise comparison out of a set of many groups may be performed by selecting the checkbox labeled ‘Compare ONLY two sample group to one another’
- 5) Review SAM options for limiting FDR cutoff, minimum and maximum genes
- 6) Optionally choose a “Name for analysis” for your results (recommended to differentiate multiple analysis runs).
- 7) Optionally add an e-mail address to be notified after the run is complete.
- 8) Click “Submit Job”

## 2.2.3 Adding data to GetExpression

In GetExpression, the data can be further analyzed, merged with different results sets and viewed in Cytoscape.

- 1) In your SAM or t-test result page, under Add Results to Get Expression → Add Data, choose **Add Data Link**
- 2) Click “Check Condition Names” to see whether this condition already exists. If the name is unique proceed to the next step, otherwise go back to the previous step and create a new name
- 3) Click the button “Upload Conditions.” Wait for the data to load and go to the GetExpression page after the data is loaded

- 4) Click the link Go to Get Expression Page **here** at the bottom of the page OR choose the **GetExpression** link on the left menu bar

## 2.3 Querying GetExpression and Launching Cytoscape

### 2.3.1 Query Data

- 1) Select the desired project(s) – Conditions from multiple projects may be chosen.
- 2) Click the **GetExpression** link on the left menu bar
- 3) Within the Conditions menu box select the conditions of interest
- 4) The Following columns **MUST** be selected in order to ensure the data can be loaded into Cytoscape:  
Data Columns to Display:
  - Log 10 Ratio
  - False Discovery RateDisplay Options:
  - Show All Conditions if one condition meets criteria
  - Pivot Conditions as columns
- 5) Click the “Query” button

### 2.3.2 Launch Cytoscape

- 6) Click the link next to Download ResultSet in Format:→ **Cytoscape**
- 7) When the web page come back, click the first three Links (The Boss must go first):
  - Boss
  - Network
  - Data Matrix Viewer







Figure 2.6 Java Web Start of Cytoscape, Boss, and Data Matrix Viewer

### 2.3.3 Viewing Expression Data in Cytoscape

- 1) Make sure all three programs are running: Gaggle Boss, Cytoscape and Data Matrix Viewer(DMV)
- 2) Load the Expression Data
  - a. Click on the DMV window
  - b. Click on the folder "Expression"
  - c. Click on the button to the right of "new" 
- 3) Load the expression data into Cytoscape
  - a. Within the DMV go up to the drop down menu currently labeled "None", and select the first Condition in the list. This should make the circle (gene nodes) change color in Cytoscape
  - b. To play a movie click the icon 
- 4) Manually Select the genes that both conditions share in common in Cytoscape
  - a. Find the Cytoscape window
  - b. Click and drag to highlight the gene nodes. Use Shift if all of the genes cannot be selected at one time
  - c. Right click to view the gene names
- 5) Use a filter to select all the genes that both conditions have in common
  - a. Click the "Grid" icon located at the top-middle of the Cytoscape menu bar
  - b. Find the tab "Topology 1"


- c. Enter 2 in the field “Minimum Number of Neighbors”
  - d. Enter 1 in the field “Within Depth”
  - e. Click Select
- 6) Graph the Selected genes in the DMV to view the global expression profile
- a. Make sure no genes are selected in the DMV. Click the “Clear” button in the DMV window
  - b. In Cytoscape with the genes of interest highlighted, click the “Broadcast” button, located in the upper right corner.
  - c. Go to the DMV window.
  - d. Select the Log10Ratio Tab. It should indicate that “7” Genes are selected.



- e. Click the icon
- 7) Find Genes that have a similar expression profile
- a. Make sure no genes are currently selected. Go to the Log10ratios window and click the “Clear” button
  - b. Within the graph made in step 5, click a gene name with an interesting profile.



- c. Click the “Broadcast” icon Located in the upper left hand part of the screen
- d. Go back to the Log10ratio screen.



- e. Click the icon
  - f. Move the slider so the “Threshold” reads 99
  - g. Graph the data by clicking the icon in 5-e
  - h. To highlight the genes in Cytoscape click the Broadcast button.
- 5-c
- i. Note that all the genes selected in the DMV may not be viewable in Cytoscape since the expression network is a small subset of all the expression data loaded from the GetExpression query.

- 8) Apply GO annotation
- a. Click the information icon top-middle of the Cytoscape menu
  - b. Click on “Go, Molecular Function, Homo sapiens”
  - c. Click on level “4”
  - d. Click “Apply Annotation to All Nodes”
  - e. Open the folder on the right half of the screen and walk through the different GO levels. Look at the Cytoscape screen to see what becomes highlighted.